

Perspective

Advancing imprinted polymers: Pre-orientation, self-assembly, and porosity in 3D printing

Amelie Huber,^{1,2} Benedikt Keitel,^{1,2} Sherman Lesly Zambou Jiokeng,² Boris Mizaikoff,^{1,2,*} and Mehmet Dinc^{1,2,*}¹Hahn-Schickard, Sedanstraße 14, 89077 Ulm, Germany²Institute of Analytical and Bioanalytical Chemistry, Ulm University, Albert-Einstein-Allee 11, 89081 Ulm, Germany*Correspondence: boris.mizaikoff@uni-ulm.de (B.M.), mehmet.dinc@hahn-schickard.de (M.D.)<https://doi.org/10.1016/j.isci.2025.114378>**SUMMARY**

Targeted isolation of specific compounds from complex matrices remains a central challenge across analytical, biomedical, and environmental sciences. Molecularly imprinted polymers (MIPs) offer a promising solution by enabling specific and reproducible binding of analytes. However, conventional imprinting techniques often suffer from heterogeneous binding sites and limited spatial control, restricting their broader utility. We see great potential in combining polymerization-induced phase-separation with photocuring 3D printing and oriented imprinting to fabricate monolithic, hierarchically porous MIPs with homogeneous and accessible binding moieties. This emerging platform enables the design of highly tailorable materials featuring distinct specificity and tunable architectures while being straightforward to produce, adapt, and scale. Beyond addressing current limitations, this strategy opens new avenues for translating imprinting technologies into industrially relevant formats. In this perspective, we outline the conceptual and technological implications of 3D-printed oriented MIPs, discuss their potential across disciplines, and highlight key challenges and opportunities for future research and application.

INTRODUCTION

In the search for a synthetic approach to mimic the specific lock-and-key principle of receptors in nature, molecularly imprinted polymers (MIPs), also known as artificial receptors, were developed.^{1–6} MIPs are highly cross-linked polymeric matrices containing tailor-made binding sites complementary to the shape, size, and distribution of functional groups of the target analyte, referred to as the template.⁷ For the fabrication of such biomimetic materials, a substance of interest is used as a structure-directing template during the copolymerization of suitable monomeric building blocks, i.e., functional monomers and cross-linkers, in a well-adjusted porogen system. After extraction of the template, the complementary binding sites remain in the polymer to yield a receptor-like structure for the template.⁸

Compared to their biological analogs, e.g., receptors or antibodies, MIPs show similar specificity toward analytes⁹ while offering higher stability in varying pH values and temperatures,⁷ being facile and usually inexpensive in preparation,¹⁰ and allowing the use of several analytes and reagents during synthesis and analyses.^{11,12} These characteristics predestine MIPs for different applications, for example, as highly specific separation substances,⁹ as recognition elements in sensors,¹³ in drug-delivery systems,¹⁴ and in catalysis.¹⁵

Molecular imprinting techniques have evolved to include not only traditional methods such as bulk polymerization and sol-gel processes but also more recent approaches like suspension, emulsion, seed, and precipitation polymerization.^{7,16–18}

MIPs can feature recognition sites within the bulk matrix, in a core-shell architecture, or confined to the surface.^{1,3,19–23}

Bulk imprinting, the classical approach, involves polymerization into monolithic structures that require grinding prior to use, with binding sites randomly distributed throughout the matrix. While this method is simple and compatible with a wide range of monomers and templates, it often suffers from slow mass transfer rates and significant template leaching, which limits specificity and overall binding capacity.^{3,24} In contrast, core-shell and surface-imprinted MIPs feature binding sites localized at or in close proximity to the surface of a supporting substrate, allowing for faster binding kinetics due to reduced diffusion distances, minimized template leaching, and enhanced binding capacities.^{1,10,21–23}

Nevertheless, these imprinting techniques often result in uncontrolled binding site formation, as the functional monomer-template complexes are fixed randomly arranged within the polymer matrix or within the polymer-shell.^{25,26} To achieve homogeneously oriented self-assembly between functional monomers and the template prior to imprinting, oriented template immobilization onto a surface is the strategy of choice (Figure 1).^{27,28}

The utilization of an orientation agent that forces the template in a specific spatial orientation prior to imprinting is referred to as oriented imprinting. After anchoring the template, the remaining functional groups are reproducibly accessible to functional monomers, enabling spatially controlled covalent or non-covalent interactions. This results in the formation of homogeneous






Figure 1. Non-oriented template immobilization in contrast to oriented template immobilization using an orientation agent

binding sites upon polymerization and ultimately improves the imprinting efficiency (Figure 2), as they offer superior binding specificity.^{27,29}

In literature, different oriented imprinting strategies using various orientation agents for template fixation are reported, such as inhibitors,³⁰ glutaraldehyde,^{31,32} boronates,^{25,27,33–38} aptamers,²⁷ and chelating agents.^{39,40} The orientation agents can be immobilized on several types of substrates, such as (nano)particles,^{28,31,32,35,38–43} 96-well microplates,³³ silicon wafers,³⁶ and monolithic capillaries.³⁶ Particle-based substrates could be based, e.g., on quantum dot (QD),³⁵ metal-organic framework (MOF),⁴² silica (SiO₂),^{32,38,39,43} or magnetic (Fe₃O₄)^{28,31,38,40,41} particles. Like their conventional counterparts, oriented MIPs have applications in various fields such as separation techniques,⁴⁴ sensing,⁴⁵ drug-delivery systems,⁴⁶ and catalysis.⁴⁷

Parallel to advances in molecular imprinting, additive manufacturing (AM), especially 3D polymer printing, gained increasing importance, driven by its exciting applications and rapid technological progress. As a key player in the manufacturing revolution, AM addresses major trends like customization and small-batch production.⁴⁸ This method allows the manufacturing of complex structures in an additive point-by-point, line-by-line, or layer-by-layer polymerization process, whereby the shape and dimensions of the object to be printed are taken from a 3D computer-aided design (CAD).^{49–52} Today, AM includes many different technologies, such as photopolymerization (also known as photocuring or photo-cross-linking) techniques like digital light processing (DLP), stereolithography (SLA), continuous liquid interface production (CLIP), and liquid crystal display (LCD) masked photocuring technology.^{50,53,54} These technologies allow fine-tuning of the physicochemical and mechanical properties of the final polymer object by selecting appropriate precursors and processing conditions. Among the notable advantages of 3D printing are its design flexibility, reduced production times, demand-driven manufacturing, and enhanced industrial sustainability through optimized processes.⁴⁸ Due to its versatile and tunable properties, 3D printing is nowadays widely used, e.g., in microfluidics,⁵⁰ biomedicine,⁵⁰ soft robotics,⁵⁵ surgery,⁵⁶ (tissue) engineering,⁵⁷ dentistry,⁵⁷ and cell cultivation.⁵⁸ The technology of 3D printing is increasingly finding its place in the field of molecular imprinting, enabling the direct fabrication of application-tailored functional materials.^{12,59–61} Among photocuring methods, LCD-

based 3D printing offers distinct advantages over point-based techniques, which, despite their high resolution, suffer from slow printing speeds and high costs. LCD-masked photocuring enables fast, layer-wide polymerization with sufficient resolution and structural fidelity. Additionally, LCD-based 3D printers are cost-effective and often easy to adapt for experimental use.⁵⁹

Recently, we fabricated 3D-MIPs by polymerization-induced phase-separation (PIPS) in a scalable and reproducible LCD-masked photocuring AM process.⁵⁹ In this study, monolithic, hierarchically porous substrates with complex macroscopic 3D designs and with excellent binding specificity toward the used template (i.e., 17 β -estradiol), were generated. The term “hierarchically porous” refers to a combination of macroscopic pores obtained by using a digital printing file and microscopic pores resulting from the use of PIPS-promoting porogens.

Broadening the use cases of oriented MIPs, we see an unrealized potential in combining this concept with PIPS in photocuring AM technology, especially LCD-masked photocuring technology.

PERSPECTIVE

Oriented 3D-printed MIPs could be fabricated by implementation of orientation agent-template composites during the polymerization of PIPS promoting resins in LCD-masked photocuring processes, leading to highly specific materials with tunable macroscopic (added 3D features in the 3D model) and inherent porosity, allowing enhanced mass transfer and spatial recognition of small and large molecules such as proteins (Figures 3 and 5).

Substrate immobilized orientation agent-template complexes can be either homogeneously dispersed in the photocurable resin containing suitable monomers and PIPS-promoting porogen(s) (as depicted in Figure 3), or directly mounted on the surface of the build platform. In both cases, self-assembly with the functional monomers occurs, forming a spatially organized pre-polymerization complex between the template and the functional monomers. Upon photocuring, the solid substrates are expected to migrate to the interfaces of the two liquid phases (polymer-rich and polymer-poor phase) according to the pickering effect, resulting in 3D-designed porous MIPs that carry the substrates at the boundaries where the pores intersect the surrounding polymer.⁶² For rebinding of the template, the composite




Figure 2. Schematic 3D representation of the production process of oriented MIPs

material must be removed, exposing the 3D-fixed binding sites adjacent to the resulting cavity complementary to the shape and size of the substrate. The 3D-printed MIP consequently could show porosity on three different levels, as it possesses macroscopic “printed” porosity (resulting from its applied CAD), a pore-network emerging from PIPS, and additional cavities with active binding sites due to removal of the orientation agent-template composite. All three types of porosity could be tailored either by adjusting the CAD, the amount and type of porogen(s), or the substrate, leading to an enhanced mass transfer and rendering this synthesis strategy even more interesting for various applications (Figure 4).

By employing this approach, direct printing of MIPs onto electrodes may become feasible, thus streamlining the process of electrochemical sensing.^{63–67} Moreover, MIPs could serve as specific target elements within microfluidic systems, enabling on-line analyte enrichment. In this scenario, the entire microfluidic setup could be printed as ready to use.^{68–70} Since the feasibility of multi-material vat photopolymerization has been established, this concept could be extended to the fabrication of MIP-assisted fully 3D-printed integrated sensors.⁵⁴ Furthermore, by performing 4D printing, thus integrating stimuli-responsiveness, it could be possible to create polymers that respond to external stimuli, such as pH, temperature, or the addition of a solvent. This could lead to physicochemical alterations in the scaffold, potentially resulting in the release or binding of a template molecule.^{71–73} Herein, we especially want to highlight the potential for epitope imprinting, as this is, up to now, a challenging field.⁷⁴ With the presented method, we expect that monolithic, hierarchically porous, oriented epitope-imprinted polymers with active binding sites that favor specific orientation for peptide binding, can be generated. These polymers can presumably rebind the corresponding, comparatively more spatially demanding biomacromolecules since the binding pockets are accessible not only through the pore network but also via the remaining voids derived upon removal of the substrate for pre-orientation (Figure 5).^{48,75–77}

A real-world example of the potential of monolithic, hierarchically porous 3D-printed oriented MIPs in the biomedical field could be the transfer to additively manufactured nasopharyngeal swabs (NPSs). During the COVID-19 pandemic local and global supply chain disruption of essential biomedical devices, such as NPSs for SARS-CoV-2 sampling, which is the reference sampling method recommended by the WHO, occurred.⁷⁸ Since 3D printing is already established in biomedical fields, the fabrication of additively manufactured biocompatible NPSs was suggested to meet the supply chain needs in addition to the mass production of conventional flock-head swabs. This technology enables rapid prototyping and cost-effective production on demand, while ensuring biocompatibility, safety, accuracy, comfort, thermal resistance, (design) flexibility, ease of handling, and functionality. Moreover, open-source collaborative design, enabled by the global exchange of digital models between academia and industry, accelerated the iterative innovation and optimization of both production and performance of 3D-printed NPSs. Studies have shown that these NPSs outperform conventional NPSs in terms of mucus collection and retention, as well as mechanical properties for greater comfort during sampling.^{78–80} We see great potential in introducing the proposed 3D printing process of oriented MIPs for virus binding into the AM of NPSs. This facilitates further functional improvements of these innovative biomedical devices by enabling affordable and accessible production of specific and selective swabs.

We anticipate that this method will have a tremendous impact on the current production of biomimetic receptors, as it could enable their simple, flexible, demand-driven, reproducible, and scalable production. The templates to be imprinted could correspond to ions, (small) molecules like, e.g., active pharmaceutical ingredients, or even epitopes for binding of larger entities, such as proteins, viruses, and bacteria, while the monomeric building blocks of the photocurable resin could be based on, e.g., acrylates or epoxides. The outstanding innovation of this approach is not only that it renders oriented imprinting easy to realize




Figure 3. Manufacturing process of porous 3D-printed MIPs with oriented binding sites

during on demand 3D printing, but in particular that it enables the production of 3D monolithic biomimetic receptors for large molecules, such as proteins or even viruses.

FUTURE PERSPECTIVES

We are confident that by combining the concepts of photocuring AM, PIPS, and oriented imprinting, hierarchically porous 3D structures with remarkable binding affinity and analyte specificity can be prepared in a straightforward, tailored, scalable, demand-driven, and reproducible manner, rendering them highly versatile for a broad range of uses. This proposed strategy demonstrates the powerful combination of mechanical performance, determined by the structure and morphology generated by AM, and the functionality, introduced by oriented MIPs in terms of analyte specificity. The innovative MIP preparation method allows new possibilities for the development of isolation processes (in both gas and liquid phase), sensors, and catalysts. Looking ahead, the integration of this methodology with emerging technologies, such as 4D printing, microfluidics, and epitope imprinting holds promise for even broader applications

in industry. Utilizing the versatility of 4D printing could enable dynamic responses to external stimuli, enhancing the functionality of MIP-based devices. Microfluidics integration may facilitate on-line specific enrichment and separation of analytes from complex matrices in a single step production method generating so-called MIP-on-a-chip devices. Furthermore, oriented epitope imprinting could be used for the specific binding of biomacromolecules or even viruses or cells on a large scale. In essence, the presented combination of these methodologies could provide groundbreaking advances in diverse fields ranging from healthcare and environmental monitoring to industrial catalysis and beyond.

Future studies should comprehensively evaluate the overall performance of 3D-printed oriented MIPs. This includes assessing and optimizing the compatibility between orientation agents and substrates across diverse templates, as well as verifying binding site homogeneity using surface or spectroscopic analyses. The mass transfer within the hierarchically porous structures should be modeled to enable rational structural design, while mechanical tests must be performed to determine the mechanical stability and robustness of the functional materials. In




Figure 4. Possible fields of application for 3D-printed oriented MIPs




Figure 5. Binding process of a protein to the corresponding binding sites within the 3D-printed hierarchically porous oriented MIP

addition, the reproducibility of the developed hierarchically porous oriented bulk MIPs throughout different batches and printing conditions should be investigated. The validation of the MIP performance should also be performed in complex real-world matrices, including for instance serum, urine, or wastewater and must be benchmarked against conventional MIP counterparts.

Application-oriented investigations should further focus on evaluating the rebinding properties and release mechanisms of target analytes under real-time sensing conditions. In biomedical applications, such as the AM of NPSs, the biocompatibility of the oriented epitope-imprinted polymers must be confirmed, analyzed, and optimized.

ACKNOWLEDGMENTS

This work was financially supported by the EU Horizon 2020 project ENVIROMED (grant no. 101057844), which focuses on developing a next-generation toolbox for greener pharmaceutical design and manufacturing to reduce environmental impact. Furthermore, it was partially funded by ZIM (Zentrales Innovationsprogramm Mittelstand) under the project MIPextract (grant KK5054610SK1).

AUTHOR CONTRIBUTIONS

Conceptualization, all authors; writing—original draft, A.H.; writing—review & editing, all authors except A.H.; visualization, A.H.; funding acquisition, B.M.; resources, B.M. and M.D.; project administration, B.M. and M.D.; supervision, B.M. and M.D.

DECLARATION OF INTERESTS

B.K., B.M., and M.D. are co-inventors on a patent application related to the LCD-based 3D-printed MIP technology described in this manuscript (PCT/EP2024/076002). One of the authors of this study, B.M., is a guest editor of the special issue “Translational and interdisciplinary advances of molecular imprinting technology” to which this paper belongs.

REFERENCES

1. Pérez, N., Whitcombe, M.J., and Vulfson, E.N. (2000). Molecularly imprinted nanoparticles prepared by core-shell emulsion polymerization. *J. Appl. Polym. Sci.* 77, 1851–1859. [https://doi.org/10.1002/1097-4628\(2000822\)77:8<1851::AID-APP23>3.0.CO;2-J](https://doi.org/10.1002/1097-4628(2000822)77:8<1851::AID-APP23>3.0.CO;2-J).

- Arshady, R., and Mosbach, K. (1981). Synthesis of substrate-selective polymers by host-guest polymerization. *Makromol. Chem.* 182, 687–692. <https://doi.org/10.1002/macp.1981.021820240>.
- Wulff, G. (1995). Molecular Imprinting in Cross-Linked Materials with the Aid of Molecular Templates— A Way towards Artificial Antibodies. *Angew. Chem. Int. Ed.* 34, 1812–1832. <https://doi.org/10.1002/anie.199518121>.
- Dickey, F.H. (1955). Specific Adsorption. *J. Phys. Chem.* 59, 695–707. <https://doi.org/10.1021/j150530a006>.
- Pauling, L. (1940). A Theory of the Structure and Process of Formation of Antibodies. *J. Am. Chem. Soc.* 62, 2643–2657. <https://doi.org/10.1021/ja01867a018>.
- Polyakov, M.V. (1931). Adsorption properties and structure of silica gel. *Zh. Fiz. Khim.* 2, 799–805.
- Lofgreen, J.E., and Ozin, G.A. (2014). Controlling morphology and porosity to improve performance of molecularly imprinted sol-gel silica. *Chem. Soc. Rev.* 43, 911–933. <https://doi.org/10.1039/C3CS60276A>.
- Zaidi, S.A. (2016). Latest trends in molecularly imprinted polymer based drug delivery systems. *RSC Adv.* 6, 88807–88819. <https://doi.org/10.1039/C6RA18911C>.
- Hu, T., Chen, R., Wang, Q., He, C., and Liu, S. (2021). Recent advances and applications of molecularly imprinted polymers in solid-phase extraction for real sample analysis. *J. Sep. Sci.* 44, 274–309. <https://doi.org/10.1002/jssc.202000832>.
- Chen, L., Xu, S., and Li, J. (2011). Recent advances in molecular imprinting technology: current status, challenges and highlighted applications. *Chem. Soc. Rev.* 40, 2922–2942. <https://doi.org/10.1039/c0cs00084a>.
- Fernandes, R.S., Dinc, M., Raimundo Jr, I.M., and Mizaikoff, B. (2017). Molecularly imprinted core-shell hybrid microspheres for the selective extraction of vanillin. *Anal. Methods* 9, 2883–2889. <https://doi.org/10.1039/C7AY00461C>.
- Paruli, E.I., Soppera, O., Haupt, K., and Gonzato, C. (2021). Photopolymerization and Photostructuring of Molecularly Imprinted Polymers. *ACS Appl. Polym. Mater.* 3, 4769–4790. <https://doi.org/10.1021/acsapm.1c00661>.
- Ozcelikay, G., Kaya, S.I., Ozkan, E., Cetinkaya, A., Nemutlu, E., Kir, S., and Ozkan, S.A. (2022). Sensor-based MIP technologies for targeted metabolomics analysis. *TrAC Trends Anal. Chem.* 146, 116487. <https://doi.org/10.1016/j.trac.2021.116487>.
- BelBruno, J.J. (2019). Molecularly Imprinted Polymers. *Chem. Rev.* 119, 94–119. <https://doi.org/10.1021/acs.chemrev.8b00171>.
- Alexander, C., Davidson, L., and Hayes, W. (2003). Imprinted polymers: artificial molecular recognition materials with applications in synthesis and catalysis. *Tetrahedron* 59, 2025–2057. [https://doi.org/10.1016/S0404-0202\(03\)00152-2](https://doi.org/10.1016/S0404-0202(03)00152-2).

16. Vasapollo, G., Sole, R.D., Mergola, L., Lazzoi, M.R., Scardino, A., Scorrano, S., and Mele, G. (2011). Molecularly Imprinted Polymers: Present and Future Prospective. *Int. J. Mol. Sci.* *12*, 5908–5945. <https://doi.org/10.3390/ijms12095908>.
17. Refaat, D., Aggour, M.G., Farghali, A.A., Mahajan, R., Wiklander, J.G., Nicholls, I.A., and Piletsky, S.A. (2019). Strategies for Molecular Imprinting and the Evolution of MIP Nanoparticles as Plastic Antibodies—Synthesis and Applications. *Int. J. Mol. Sci.* *20*, 6304. <https://doi.org/10.3390/ijms20246304>.
18. Sajini, T., and Mathew, B. (2021). A brief overview of molecularly imprinted polymers: Highlighting computational design, nano and photo-responsive imprinting. *Talanta Open* *4*, 100072. <https://doi.org/10.1016/j.talo.2021.100072>.
19. Katz, A., and Davis, M.E. (2000). Molecular imprinting of bulk, microporous silica. *Nature* *403*, 286–289. <https://doi.org/10.1038/35002032>.
20. Pérez-Moral, N., and Mayes, A.G. (2004). Noncovalent Imprinting in the Shell of Core–Shell Nanoparticles. *Langmuir* *20*, 3775–3779. <https://doi.org/10.1021/la0356755>.
21. Dong, C., Shi, H., Han, Y., Yang, Y., Wang, R., and Men, J. (2021). Molecularly imprinted polymers by the surface imprinting technique. *Eur. Polym. J.* *145*, 110231. <https://doi.org/10.1016/j.eurpolymj.2020.110231>.
22. Niu, M., Pham-Huy, C., and He, H. (2016). Core-shell nanoparticles coated with molecularly imprinted polymers: a review. *Microchim. Acta* *183*, 2677–2695. <https://doi.org/10.1007/s00604-016-1930-4>.
23. Wan, L., Chen, Z., Huang, C., and Shen, X. (2017). Core-shell molecularly imprinted particles. *TrAC Trends Anal. Chem.* *95*, 110–121. <https://doi.org/10.1016/j.trac.2017.08.010>.
24. Alexander, C., Andersson, H.S., Andersson, L.I., Ansell, R.J., Kirsch, N., Nicholls, I.A., O'Mahony, J., and Whitcombe, M.J. (2006). Molecular imprinting science and technology: a survey of the literature for the years up to and including 2003. *J. Mol. Recognit.* *19*, 106–180. <https://doi.org/10.1002/jmr.760>.
25. Xing, R., Wang, S., Bie, Z., He, H., and Liu, Z. (2017). Preparation of molecularly imprinted polymers specific to glycoproteins, glycans and monosaccharides via boronate affinity controllable-oriented surface imprinting. *Nat. Protoc.* *12*, 964–987. <https://doi.org/10.1038/nprot.2017.015>.
26. Cabaleiro-Lago, C., Hasterok, S., Gjörlöf Wingren, A., and Tassidis, H. (2023). Recent Advances in Molecularly Imprinted Polymers and Their Disease-Related Applications. *Polymers* *15*, 4199. <https://doi.org/10.3390/polym15214199>.
27. Kalecki, J., Iskierko, Z., Cieplak, M., and Sharma, P.S. (2020). Oriented Immobilization of Protein Templates: A New Trend in Surface Imprinting. *ACS Sens.* *5*, 3710–3720. <https://doi.org/10.1021/acssensors.0c01634>.
28. Dinc, M., Esen, C., and Mizaikoff, B. (2019). Recent advances on core-shell magnetic molecularly imprinted polymers for biomacromolecules. *TrAC Trends Anal. Chem.* *114*, 202–217. <https://doi.org/10.1016/j.trac.2019.03.008>.
29. He, S., Zhang, L., Bai, S., Yang, H., Cui, Z., Zhang, X., and Li, Y. (2021). Advances of molecularly imprinted polymers (MIP) and the application in drug delivery. *Eur. Polym. J.* *143*, 110179. <https://doi.org/10.1016/j.eurpolymj.2020.110179>.
30. Dinc, M., Basan, H., Diemant, T., Behm, R.J., Lindén, M., and Mizaikoff, B. (2016). Inhibitor-assisted synthesis of silica-core microbeads with pepsin-imprinted nanoshells. *J. Mater. Chem. B* *4*, 4462–4469. <https://doi.org/10.1039/C6TB00147E>.
31. Sullivan, M.V., Stockburn, W.J., Hawes, P.C., Mercer, T., and Reddy, S.M. (2021). Green synthesis as a simple and rapid route to protein modified magnetic nanoparticles for use in the development of a fluorometric molecularly imprinted polymer-based assay for detection of myoglobin. *Nanotechnology* *32*, 095502. <https://doi.org/10.1088/1361-6528/abc2d>.
32. Zhang, S., Liu, Z., Jin, S., Bai, Y., Feng, X., and Fu, G. (2021). A method for synthesis of oriented epitope-imprinted open-mouthed polymer nanocapsules and their use for fluorescent sensing of target protein. *Talanta* *234*, 122690. <https://doi.org/10.1016/j.talanta.2021.122690>.
33. Bi, X., and Liu, Z. (2014). Facile Preparation of Glycoprotein-Imprinted 96-Well Microplates for Enzyme-Linked Immunosorbent Assay by Boronate Affinity-Based Oriented Surface Imprinting. *Anal. Chem.* *86*, 959–966. <https://doi.org/10.1021/ac403736y>.
34. Bie, Z., and Chen, Y. (2021). Selective analysis of interferon-alpha in human serum with boronate affinity oriented imprinting based plastic antibody. *Talanta* *230*, 122338. <https://doi.org/10.1016/j.talanta.2021.122338>.
35. Li, P., and Liu, Z. (2024). Glycan-specific molecularly imprinted polymers towards cancer diagnostics: merits, applications, and future perspectives. *Chem. Soc. Rev.* *53*, 1870–1891. <https://doi.org/10.1039/D3CS00842H>.
36. Wang, S., Ye, J., Bie, Z., and Liu, Z. (2014). Affinity-tunable specific recognition of glycoproteins via boronate affinity-based controllable oriented surface imprinting. *Chem. Sci.* *5*, 1135. <https://doi.org/10.1039/c3sc52986j>.
37. Xu, J., Miao, H., Wang, J., and Pan, G. (2020). Molecularly Imprinted Synthetic Antibodies: From Chemical Design to Biomedical Applications. *Small* *16*, 1906644. <https://doi.org/10.1002/sml.201906644>.
38. Zhou, J., Mujahid Ali, M., Yu, W., Cheng, X., Gao, Y., and Hu, L. (2024). Oriented docking of the template for improved imprinting efficiency toward peptide with modifications. *Anal. Chim. Acta* *1301*, 342450. <https://doi.org/10.1016/j.aca.2024.342450>.
39. Xu, J., Miao, H., Zou, L., Tse Sum Bui, B., Haupt, K., and Pan, G. (2021). Evolution of Molecularly Imprinted Enzyme Inhibitors: From Simple Activity Inhibition to Pathological Cell Regulation. *Angew. Chem.* *133*, 24731–24738. <https://doi.org/10.1002/ange.202106657>.
40. Zhu, R., Wang, X., Ma, J., and Jia, Q. (2024). Preparation of poly(caffeic acid)-coated epitope molecularly imprinted polymers and investigation of adsorption performance toward ovalbumin. *J. Chromatogr. A* *1716*, 464635. <https://doi.org/10.1016/j.chroma.2024.464635>.
41. Hou, H., Jin, Y., Xu, K., Sheng, L., Huang, Y., and Zhao, R. (2021). Selective recognition of a cyclic peptide hormone in human plasma by hydrazone bond-oriented surface imprinted nanoparticles. *Anal. Chim. Acta* *1154*, 338301. <https://doi.org/10.1016/j.aca.2021.338301>.
42. Hu, W., Pei, F., Du, B., Wang, J., Liang, M., Yang, L., Liu, B., Mu, X., and Tong, Z. (2024). A fluorescence-electrochemistry dual-mode imprinted sensing platform constructed by boric acid-functionalized MOF@COF core-shell composite for sensitive detection of glycoprotein. *Sens. Actuators B Chem.* *407*, 135494. <https://doi.org/10.1016/j.snb.2024.135494>.
43. Nishino, H., Huang, C.S., and Shea, K.J. (2006). Selective Protein Capture by Epitope Imprinting. *Angew. Chem. Int. Ed.* *45*, 2392–2396. <https://doi.org/10.1002/anie.200503760>.
44. Lai, S., Ouyang, X., Cai, C., Xu, W., Chen, C., and Chen, X. (2017). Surface-imprinted microspheres prepared by a template-oriented method for the chiral separation of amlodipine. *J. Sep. Sci.* *40*, 1869–1876. <https://doi.org/10.1002/jssc.201700076>.
45. Jolly, P., Tamboli, V., Harniman, R.L., Estrela, P., Allender, C.J., and Bowen, J.L. (2016). Aptamer-MIP hybrid receptor for highly sensitive electrochemical detection of prostate specific antigen. *Biosens. Bioelectron.* *75*, 188–195. <https://doi.org/10.1016/j.bios.2015.08.043>.
46. Lawai, V., Ngaini, Z., Farooq, S., Wahi, R., and Bhawani, S.A. (2024). Current advances in molecularly imprinted polymers and their release mechanisms in drug delivery systems. *Polym. Adv. Technol.* *35*, e6317. <https://doi.org/10.1002/pat.6317>.
47. Zhang, G., Yu, Y., Zhang, L., Lin, B., Wang, Y., Guo, M., and Cao, Y. (2020). Precise detection of prostate specific antigen in serum: A surface molecularly imprinted sensor based on novel cooperated signal amplification strategy. *Sens. Actuators B Chem.* *302*, 126998. <https://doi.org/10.1016/j.snb.2019.126998>.

48. Keitel, B., Dietl, S., Philipp, T., Neusser, G., Kranz, C., Sobek, H., Mizaikoff, B., and Dinc, M. (2025). 3D-Printed Highly Porous Functional Materials for the Efficient Removal of Adenovirus. *Adv. Mater. Technol.* *10*, 2401178. <https://doi.org/10.1002/admt.202401178>.
49. Rajpurohit, S.R., and Dave, H.K. (2020). Tensile Strength of 3D Printed PLA Part. In *Advances in Additive Manufacturing and Joining Lecture Notes on Multidisciplinary Industrial Engineering*, M.S. Shunmugam and M. Kanthababu, eds. (Springer Singapore), pp. 103–114. https://doi.org/10.1007/978-981-32-9433-2_8.
50. Bagheri, A., and Jin, J. (2019). Photopolymerization in 3D Printing. *ACS Appl. Polym. Mater.* *1*, 593–611. <https://doi.org/10.1021/acsapm.8b00165>.
51. Peerzada, M., Abbasi, S., Lau, K.T., and Hameed, N. (2020). Additive Manufacturing of Epoxy Resins: Materials, Methods, and Latest Trends. *Ind. Eng. Chem. Res.* *59*, 6375–6390. <https://doi.org/10.1021/acs.iecr.9b06870>.
52. Wendel, B., Rietzel, D., Kühnlein, F., Feulner, R., Hülder, G., and Schmachtenberg, E. (2008). Additive Processing of Polymers. *Macromol. Mater. Eng.* *293*, 799–809. <https://doi.org/10.1002/mame.200800121>.
53. Hu, Y., Luo, Z., and Bao, Y. (2025). Trends in Photopolymerization 3D Printing for Advanced Drug Delivery Applications. *Biomacromolecules* *26*, 85–117. <https://doi.org/10.1021/acs.biomac.4c01004>.
54. Subedi, S., Liu, S., Wang, W., Naser Shovon, S.M.A., Chen, X., and Ware, H.O.T. (2024). Multi-material vat photopolymerization 3D printing: a review of mechanisms and applications. *NPJ Adv. Manuf.* *1*, 9. <https://doi.org/10.1038/s44334-024-00005-w>.
55. Wallin, T.J., Pikul, J., and Shepherd, R.F. (2018). 3D printing of soft robotic systems. *Nat. Rev. Mater.* *3*, 84–100. <https://doi.org/10.1038/s41578-018-0002-2>.
56. Pugliese, L., Marconi, S., Negrello, E., Mauri, V., Peri, A., Gallo, V., Auricchio, F., and Pietrabissa, A. (2018). The clinical use of 3D printing in surgery. *Updates Surg.* *70*, 381–388. <https://doi.org/10.1007/s13304-018-0586-5>.
57. Quan, H., Zhang, T., Xu, H., Luo, S., Nie, J., and Zhu, X. (2020). Photocuring 3D printing technique and its challenges. *Bioact. Mater.* *5*, 110–115. <https://doi.org/10.1016/j.bioactmat.2019.12.003>.
58. Schmidleithner, C., Malferrari, S., Palgrave, R., Bomze, D., Schwentenwein, M., and Kalaskar, D.M. (2019). Application of high resolution DLP stereolithography for fabrication of tricalcium phosphate scaffolds for bone regeneration. *Biomed. Mater.* *14*, 045018. <https://doi.org/10.1088/1748-605X/ab279d>.
59. Keitel, B., Batista, A.D., Schimana, S., Mizaikoff, B., and Dinc, M. (2024). Emulsion-Free 3D Printing of Inherently Porous Molecularly Imprinted Polymers with Tailored Macroscopic Geometries. *ACS Appl. Polym. Mater.* *6*, 3690–3695. <https://doi.org/10.1021/acsapm.3c02744>.
60. Gomez, L.P.C., Spangenberg, A., Ton, X.A., Fuchs, Y., Bokeloh, F., Malval, J.P., Tse Sum Bui, B., Thuau, D., Ayela, C., Haupt, K., and Soppera, O. (2016). Rapid Prototyping of Chemical Microsensors Based on Molecularly Imprinted Polymers Synthesized by Two-Photon Stereolithography. *Adv. Mater.* *28*, 5931–5937. <https://doi.org/10.1002/adma.201600218>.
61. Conrad, P.G., Nishimura, P.T., Aherne, D., Schwartz, B.J., Wu, D., Fang, N., Zhang, X., Roberts, M.J., and Shea, K.J. (2003). Functional Molecularly Imprinted Polymer Microstructures Fabricated Using Microstereolithography. *Adv. Mater.* *15*, 1541–1544. <https://doi.org/10.1002/adma.200304602>.
62. Dubinsky, S., Petukhova, A., Gourevich, I., and Kumacheva, E. (2010). Hybrid porous material produced by polymerization-induced phase separation. *Chem. Commun.* *46*, 2578–2580. <https://doi.org/10.1039/b924373a>.
63. Khosrokhavar, R., Motaharian, A., Milani Hosseini, M.R., and Mohammad-sadegh, S. (2020). Screen-printed carbon electrode (SPCE) modified by molecularly imprinted polymer (MIP) nanoparticles and graphene nanosheets for determination of sertraline antidepressant drug. *Microchem. J.* *159*, 105348. <https://doi.org/10.1016/j.microc.2020.105348>.
64. McClements, J., Seumo Tchekwagep, P.M., Vilela Strapazon, A.L., Canfarotta, F., Thomson, A., Czulak, J., Johnson, R.E., Novakovic, K., Losada-Pérez, P., Zaman, A., et al. (2021). Immobilization of Molecularly Imprinted Polymer Nanoparticles onto Surfaces Using Different Strategies: Evaluating the Influence of the Functionalized Interface on the Performance of a Thermal Assay for the Detection of the Cardiac Biomarker Troponin I. *ACS Appl. Mater. Interfaces* *13*, 27868–27879. <https://doi.org/10.1021/acsami.1c05566>.
65. Mostafiz, B., Bigdeli, S.A., Banan, K., Afsharara, H., Hatamabadi, D., Mousavi, P., Hussain, C.M., Keçili, R., and Ghorbani-Bidkorbeh, F. (2021). Molecularly imprinted polymer-carbon paste electrode (MIP-CPE)-based sensors for the sensitive detection of organic and inorganic environmental pollutants: A review. *Trends Environ. Anal. Chem.* *32*, e00144. <https://doi.org/10.1016/j.teac.2021.e00144>.
66. Chu, T., Park, S., and Fu, K.K. (2021). 3D printing-enabled advanced electrode architecture design. *Carbon Energy* *3*, 424–439. <https://doi.org/10.1002/cey2.114>.
67. Ahmed, A.H.M.T., Naskar, H., Banerjee, S., Ghatak, B., Das, N., Tudu, B., and Bandyopadhyay, R. (2022). Electrochemical sensor based on molecularly imprinted polymer embedded graphite electrode for detecting curcumin. *Sens. Actuators Phys.* *344*, 113748. <https://doi.org/10.1016/j.sna.2022.113748>.
68. Azhdary, P., Janfaza, S., Fardindoost, S., Tasnim, N., and Hoorfar, M. (2023). Highly selective molecularly imprinted polymer nanoparticles (MIP NPs)-based microfluidic gas sensor for tetrahydrocannabinol (THC) detection. *Anal. Chim. Acta* *1278*, 341749. <https://doi.org/10.1016/j.aca.2023.341749>.
69. Burnage, S.C., Bell, J., Wan, W., Kislenko, E., and Rurack, K. (2023). Combining a hybrid chip and tube microfluidic system with fluorescent molecularly imprinted polymer (MIP) core-shell particles for the derivatization, extraction, and detection of peptides with N-terminating phosphorylated tyrosine. *Lab Chip* *23*, 466–474. <https://doi.org/10.1039/D2LC00955B>.
70. Uka, B., Kieninger, J., Urban, G.A., and Weltin, A. (2021). Electrochemical Microsensor for Microfluidic Glyphosate Monitoring in Water Using MIP-Based Concentrators. *ACS Sens.* *6*, 2738–2746. <https://doi.org/10.1021/acssensors.1c00884>.
71. Andreu, A., Su, P.-C., Kim, J.-H., Ng, C.S., Kim, S., Kim, I., Lee, J., Noh, J., Subramanian, A.S., and Yoon, Y.-J. (2021). 4D printing materials for vat photopolymerization. *Addit. Manuf.* *44*, 102024. <https://doi.org/10.1016/j.addma.2021.102024>.
72. Kocak, G., Tuncer, C., and Bütün, V. (2017). pH-Responsive polymers. *Polym. Chem.* *8*, 144–176. <https://doi.org/10.1039/C6PY01872F>.
73. Kitayama, Y., and Isomura, M. (2020). Molecularly imprinted polymer particles with gas-stimuli responsive affinity toward target proteins prepared using switchable functional monomer. *Polymer* *203*, 122781. <https://doi.org/10.1016/j.polymer.2020.122781>.
74. Dietl, S., Sobek, H., and Mizaikoff, B. (2021). Epitope-imprinted polymers for biomacromolecules: Recent strategies, future challenges and selected applications. *TrAC Trends Anal. Chem.* *143*, 116414. <https://doi.org/10.1016/j.trac.2021.116414>.
75. Amorim, M.S., Sales, M.G.F., and Frasco, M.F. (2022). Recent advances in virus imprinted polymers. *Biosens. Bioelectron. X* *10*, 100131. <https://doi.org/10.1016/j.biosx.2022.100131>.
76. Yang, K., Li, S., Liu, L., Chen, Y., Zhou, W., Pei, J., Liang, Z., Zhang, L., and Zhang, Y. (2019). Epitope Imprinting Technology: Progress, Applications, and Perspectives toward Artificial Antibodies. *Adv. Mater.* *31*, 1902048. <https://doi.org/10.1002/adma.201902048>.
77. Zhang, W., Zhang, Y., Wang, R., Zhang, P., Zhang, Y., Randell, E., Zhang, M., and Jia, Q. (2022). A review: Development and application of surface molecularly imprinted polymers toward amino acids, peptides, and proteins. *Anal. Chim. Acta* *1234*, 340319. <https://doi.org/10.1016/j.aca.2022.340319>.

78. Van Der Elst, L.A., Gokce Kurtoglu, M., Leffel, T., Zheng, M., and Gumenik, A. (2020). Rapid Fabrication of Sterile Medical Nasopharyngeal Swabs by Stereolithography for Widespread Testing in a Pandemic. *Adv. Eng. Mater.* 22, 2000759. <https://doi.org/10.1002/adem.202000759>.
79. Arjunan, A., Zahid, S., Baroutaji, A., and Robinson, J. (2021). 3D printed auxetic nasopharyngeal swabs for COVID-19 sample collection. *J. Mech. Behav. Biomed. Mater.* 114, 104175. <https://doi.org/10.1016/j.jmbbm.2020.104175>.
80. Xiao, R., Li, J., Wong, Y.F., Chu, L.T., Zhu, Y., Wang, L., Wu, M., Zhang, D., Gong, M., Lai, J., et al. (2024). Enhancing detection accuracy via controlled release of 3D-printed microlattice nasopharyngeal swabs. *Commun. Eng.* 3, 40. <https://doi.org/10.1038/s44172-024-00185-5>.